

PR DETECT™

Dengue Rapid Test Series



Dengue Rapid Test Series, is an *in vitro* qualitative immunochromatographic assay for the rapid detection of dengue NS1 antigen and anti-dengue IgG/IgM antibodies in human serum, plasma or whole blood specimens.



FAST
Test result in 10-20 mins.



FLEXIBLE
Available in 3 formats - NS1 Ag, IgG/IgM, NS1 Ag & IgG/IgM.






CONVENIENT
Rapid, 24 months shelf life at room temperature.



HIGHLY INFORMATIVE
NS1 & IgG/IgM - Detection of dengue NS1 antigen and IgG/IgM antibodies of all serotypes (DEN 1-4).

PR DETECT™	Combined Diagnostics for Dengue NS1 Antigen and IgG/IgM Antibodies		
	Sensitivity	Specificity	Specimen
NS1 Ag	95.8%	96.1%	Whole blood/ Serum/ Plasma
IgG/IgM	95.7%	99.9%	Whole blood/ Serum/ Plasma

Clinical Stages




ORDERING INFORMATION

Catalog No.	Product Name	Format	Width	Specimen	Tests/ Kit
PR-DEN-1	ProDetect™ Dengue NS1 Ag Rapid Test*	Cassette	4mm	Whole blood/ Serum/ Plasma	25
PR-DEN-2	ProDetect™ Dengue IgG/IgM Rapid Test*	Cassette	4mm	Whole blood/ Serum/ Plasma	25
PR-DEN-3-1	ProDetect™ Dengue Duo NS1 Ag & IgG/IgM Rapid Test*	Cassette	4mm	Whole blood/ Serum/ Plasma	10
PR-DEN-3-2	ProDetect™ Dengue Duo NS1 Ag & IgG/IgM Rapid Test*	Cassette	4mm	Whole blood/ Serum/ Plasma	25

*Products are CE-IVD certified.



For Dengue NS1 Antigen Test



For Dengue IgG/IgM Antibodies Test

