



Multi-drug Screen Test Cup

A rapid, one step screening test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human urine.

This multi-drug screen test cup is a competitive immunoassay utilizing highly specific reactions between antibodies and antigens for the detection of multiple drugs and drug metabolites in human urine without the use of an instrument.



FAST

Test result in 5-10 mins



FLEXIBLE

Different combination of drug panels available.




CONVENIENT

Rapid, stable at room temperature.



HIGHLY INFORMATIVE

Temperature strips are attached to the vial to help read result of the test.



ORDERING INFORMATION

| Catalog No. | Product Name | Combination | Cups/ Box |
|-------------|---|--------------------------|-----------|
| PR-KC-6 | ProDetect™ Push Button 6 in 1 Multi-Drug Test Cup | AMP/COC/OPI/BZO/THC/MDMA | 25 |