

# PR DETECT™ Zika Rapid Test



Zika virus spreads to people primarily through the bite of an infected *Aedes* species mosquito (*Aedes aegypti* and *Aedes albopictus*). Zika can also be passed through sex from a person who has Zika to his or her sex partners and it can be spread from a pregnant woman to her fetus.

Zika outbreaks are currently happening in many countries and territories. The mosquitoes can become infected with Zika virus and spread the infection to many parts of the world.

Use **ProDetect™ Zika Rapid Test** for rapid and convenient immunochromatographic *in vitro* assay for the differential detection of IgG and IgM antibodies to Zika virus in human serum, plasma and/or whole blood samples.



## FAST

Test result in 10-20 mins.



## FLEXIBLE

Qualitative detection of Zika virus infection.



## CONVENIENT

Rapid, 24 months shelf life at room temperature.



## HIGHLY INFORMATIVE

Enables early screening and rapid diagnosis of Zika virus infection.

## Countries with active Zika transmission





Barbados	Honduras
Bolivia	Martinique
Brazil	Mexico
Colombia	Panama
Ecuador	Paraguay
El Salvador	Puerto Rico
French Guiana	Saint Martin
Guadeloupe	Suriname
Guatemala	Venezuela
Guyana	Samoa
Haiti	Cape Verde

Centers for Disease Control and Prevention

### Zika Symptoms

Many people infected with Zika won't have symptoms or will only have mild symptoms. The most common symptoms are fever, rash, joint pain, or red eyes. Other common symptoms include muscle pain and headache. Symptoms can last for several days to a week. People usually don't get sick enough to go to the hospital, and they very rarely die of Zika. Once a person has been infected with Zika, they are likely to be protected from future infections.

Source: [www.cdc.gov/zika](http://www.cdc.gov/zika)



## ORDERING INFORMATION

Catalog No.	Product Name	Format	Specimen	Tests/ Kit
PR-ZIKA	ProDetect™ Zika Rapid Test*	Cassette	Whole blood/ Serum/ Plasma	25

\*Product is CE-IVD certified.